

AS120

GC & GC/MS AUTOSAMPLER

A reliable and affordable solution



KEY FEATURES

- Direct syringe injection
- Dual inlet injection capability
- Variable injection and syringe fill rates
- Method linking
- Multiple solvent/waste bottles
- Syringe size 1.5 µl-100µl
- Multiple samples per vial
- Air gapping capability
- Compatible with any GC
- Multiple tray options:
 - 120 position 2ml vials
 - 60 position 10 ml vials
- Uses off-the-shelf syringes

EST 
analytical *Rely On Us* 



SPECIFICATIONS

Size: Approx 25" L X 12" W X 17" H

Weight: Approx 17.3 lbs.

Tray Options:

- 1.8 ml X 120 position
- 10 ml x 60 position

Power In: 100-240 VAC; 1.5A 50-60 Hz

OVERALL VALUE

- Best Value on the Market
- Proudly made in the USA

BUY WITH CONFIDENCE

- EST analytical is a worldwide leader in GC sample introduction

EST  **Rely On Us**
analytical 